

## ScienceDesk™ Accessories



U.S. Patent: D455,300S  
GB Patent: GB2.365.991B

### Accessories Include:

- Shelves
- Keyboard Holders
- Wrist Rest
- Flat Screen Support Brackets
- Post Assemblies
- Earth Grounding Kits

The accessories for Thorlabs' ScienceDesk series of workstations provide an array of custom configurations to meet almost all workspace requirements. A multitude of shelves are available, including overhead shelves, side shelves, and lower shelves to support test instrumentation, computers, and tools. All shelves are made using ESD materials. A monitor mount allows standard flat panel monitors to be supported off the work surface. Keyboard mounts and a wrist rest are also available to enhance an ergonomic workspace further. Earth grounding kits make the ScienceDesk workstations a great choice for more sensitive applications.

Our ScienceDesk accessories are available for all models of the ScienceDesk frames. Accessories can be mounted in a multitude of locations to maximize your comfort and productivity, thus allowing a workspace to be reconfigured as requirements change.

## Under Shelf



The under shelf spans the rear of the ScienceDesk, providing a 260mm (10.2") shelf depth. The shelf is mounted on the side rails of the lower frame.

The top surface of the shelf is coated with electrostatic-resistant material to minimize electrostatic discharge, and it has a 75kg (165lbs) load capacity.

### Under Shelf: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK LENGTH	SHELF DIMENSION
PSX300	\$ 161.50	£101.70	€ 150.20	¥ 1,542.30	750mm	910 x 260mm (33.5" x 10.2")
PSX301	\$ 187.00	£117.80	€ 173.90	¥ 1,785.90	900mm	1060 x 260mm (41.7" x 10.2")
PSX302	\$ 224.00	£141.10	€ 208.30	¥ 2,139.20	1200mm	1360 x 260mm (53.5" x 10.2")