

Tables & Breadboards

Overhead Shelf With Monitor Mount



The rear-mount overhead shelf with monitor mounting brackets spans the rear of the ScienceDesk™, providing a shelf width of 300mm (11.8"). The shelf is mounted on two 600mm posts (included), and the shelf height is adjustable. The top surface is made with ESD materials to minimize electrostatic discharge. Each shelf has a 75kg (165lbs) load capacity and comes with two mounting brackets with a multidirectional pivot clamp, which supports most types of flat screens.

The side mount shelf is designed to span the length of the ScienceDesk as opposed to the rear mount shelf (shown left), which spans the width of the ScienceDesk.

Rear Mount Overhead Shelves: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK SIZE	SHELF DIMENSION
PSX210	\$ 770.10	£ 485.20	€ 716.20	¥ 7,354.50	750mm	850mm x 300mm (33.5" x 11.8")
PSX211	\$ 785.40	£ 494.80	€ 730.40	¥ 7,500.60	900mm	1000mm x 300mm (39.4" x 11.8")
PSX212	\$ 856.80	£ 539.80	€ 796.80	¥ 8,182.40	1200mm	1300mm x 300mm (51.2" x 11.8")

Side Mount Overhead Shelves: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK SIZE	SHELF DIMENSION
PSX213	\$ 732.70	£ 461.60	€ 681.40	¥ 6,997.30	750mm	600mm x 300mm (33.5" x 11.8")
PSX214	\$ 753.10	£ 474.50	€ 700.40	¥ 7,192.10	900mm	800mm x 300mm (39.4" x 11.8")

Please contact your local sales office for a quotation including shipping.

ScienceDesk™ Overhead Shelves

The standard overhead shelf spans the rear of the ScienceDesk™, providing a shelf with a width of 300mm (11.8"). Two models of the overhead shelf can be mounted to either side of the ScienceDesk frame. The shelf is mounted on two 600mm posts (included) and the height is adjustable.

- Shelf Load Capacity is 75kg (165lbs)
- Made With ESD Materials
- Two Mounting Posts Included

Standard Overhead Shelves: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK LENGTH	SHELF DIMENSION
PSX200	\$ 544.00	£ 342.70	€ 505.90	¥ 5,195.20	750mm	850mm x 300mm (33.5" x 11.8")
PSX201	\$ 569.50	£ 358.80	€ 529.60	¥ 5,438.70	900mm	1000mm x 300mm (39.4" x 11.8")
PSX202	\$ 632.40	£ 398.40	€ 588.10	¥ 6,039.40	1200mm	1300mm x 300mm (51.2" x 11.8")

Side Overhead Shelf: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK LENGTH	SHELF DIMENSION
PSX203	\$ 527.00	£ 332.00	€ 490.10	¥ 5,032.90	750mm	850mm x 300mm (33.5" x 11.8")
PSX204	\$ 544.00	£ 342.70	€ 505.90	¥ 5,195.20	900mm	1000mm x 300mm (39.4" x 11.8")

Single Shelf

The single shelf can be mounted on the upper side or rear rails of the ScienceDesk base. It provides a shelf size of 300 x 300mm (11.8" x 11.8"). The top surface is coated with electrostatic-dissipative material to minimize electrostatic discharge.

- The Single Shelf has a 30kg (66lbs) Load Capacity
- Made With ESD Materials
- One 600mm Mounting Post Included

Single Shelf with Post: In Stock

ITEM#	\$	£	€	RMB	SHELF DIMENSION
PSX190	\$ 340.00	£ 214.20	€ 316.20	¥ 3,247.00	300mm x 300mm (11.8" x 11.8")



Visit
www.thorlabs.com
For Mechanical
Drawings and
Our New
Solid Models

