Motion Control Motorized Actuators – DC Servo

Single Axis Stages Multi-Axis Stages

Flexure Stage Accessories Motorized

Rotation Stages

Drive Electronics & Auto-Alignment

Actuators & Adjusters

Brief Tutorials





Our Z7 series motorized micrometers are engineered for use with optical positioning devices. The high-resolution design in a lightweight and compact package makes these actuators ideally suited for demanding optical laboratory automation applications. Additional features include an optical encoder for closed-loop operation as well as built-in limit switches that provide overdrive protection.

The Z706 model has been designed specifically for use in mirror mounts that have 1/4"-80 manual adjusters. Just remove the existing 1/4"-80 drive screw and thread in our Z706 actuator. See page 352 for details of our T-Cube motor driver.

- Motorized Actuator, 12V DC Servo Motor
- Mirror Mount Drop-In Replacement for 1/4"-80 Actuators
- Actuator Travel: 6mm, 12mm, and 25mm
- 0.05µm Minimum Incremental Movement
- Backlash <8μm (Before Software Compensation)
- Compatible With Many Standard Controllers/Drivers
- Limit Switches Prevent Accidental Actuator Damage
- High Precision Rotary Encoder (48 Count Encoder x 256:1 Gear Head = 12,288 Counts per Revolution)

Z706

6mm Travel (1/4"-80 Barrel)



12mm Vacuum Rated Actuator, 3/8" Barrel



¥4,679.50

Z725B

Z612BV \$490.00

25mm Travel (3/8" Diameter Barrel)

£308.70

€455,70

ITEM#	\$	£	€	RMB	DESCRIPTION
Z725B	\$ 490.00	£ 308.70	€ 455,70	¥ 4,679.50	25mm Motorized Actuator, 3/8" Barrel
Z625BV	\$ 500.00	£ 315.00	€ 465,00	¥ 4,775.00	25mm Vacuum Rated Actuator, 3/8"Barrel

Vacuum Option

Our new vacuum compatible, motorized actuators offer the same great functionality as our standard actuators but with vacuum rated servo-motors, phosphorus bronze internal coupling mechanism, and high vacuum grease (see price boxes).

- Vacuum Rating to 10⁻⁶ Torr
- Ø3/8" or 1/4"-80 Mounting Barrel
- See Web for Details



THORLADS