


AFM1 - May 14, 2015

Item # AFM1 was discontinued on May 14, 2015. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

LAB MATS

- ▶ **Multilayered Adhesive Mat with Numbered Layers (30 Total)**
- ▶ **Interlocking Grease- and Chemical-Resistant Rubber Mats**



AFM1
Two Mats Being
Joined Together

ESD2436
Multilayered Adhesive Mats




[Hide Overview](#)

OVERVIEW

Thorlabs offers two solutions for floor mats in the laboratory: anti-fatigue and clean-walk mats. The black, rubber, anti-fatigue mats can be fit together to form various lengths and offer a chemical- and grease-resistant floor solution. Our Clear-Walk adhesive floor mats have a tacky adhesive-coated film that effectively captures dirt and dust from foot traffic and equipment wheels before entering controlled environments.

[Hide Anti-Fatigue Mats](#)

Anti-Fatigue Mats



These black, rubber, anti-fatigue mats resist oil, grease, chemicals, and acids. Their interlocking design allows for end sections to be joined together or for two end sections to be combined with several center sections to form any length in 900 mm (35.43") increments. A single mat will span the width of a 3' (900 mm) ScienceDesk, while two mats joined together will span nearly the length of a 6' (2 m) optical table.

A waffle-shaped backing cushions the feet and absorbs shock, while a textured surface provides traction and prevents slipping. All exposed edges are tapered to minimize trip hazards and to allow trolleys and carts to roll over without snagging. These mats are not conductive and

may be subject to static electricity build up.

Part Number	Description	Price	Availability
AFM1	Anti-Fatigue Mat, 900 mm x 700 mm (35.43" x 27.56"), End Piece	\$37.00	2 Weeks
AFM2	Anti-Fatigue Mat, 900 mm x 700 mm (35.43" x 27.56"), Center Piece	\$37.00	2 Weeks
AFM3	Anti-Fatigue Mat, 900 mm x 700 mm (35.43" x 27.56"), Individual Piece	\$37.00	2 Weeks

Hide Clean-Walk Mats

Clean-Walk Mats



These adhesive floor mats have a tacky adhesive-coated film that effectively captures dirt and dust from foot traffic and equipment wheels before entering controlled environments. Each layer is numbered, which helps both to ensure one-sheet removal and to identify the number of remaining sheets. Low-profile construction reduces the risk of tripping and allows free movement of foot or wheel traffic.

The adhesive mat contains 30 transparent sheets, measuring 24" x 36" (609.6 mm x 914.4 mm), with polyethylene film backing. Simply remove the clear backing, marked with a piece of tape, from the back of the mat to expose a tacky adhesive backing that sticks directly to a

clean, dry floor.

Part Number	Description	Price	Availability
ESD2436	24" x 36" (610 mm x 910 mm) Adhesive Mat, 30 Sheets	\$52.10	Today

Visit the *Lab Mats* page for pricing and availability information:

https://www.thorlabs.com/newgrouppage9.cfm?objectgroup_id=2862