



Translation Stages

Motorized Mirror Mounts

Rotation Stages

3 Axis Flexure Stages & Accessories

4, 5 & 6-Axis Stages

Drive Electronics & Auto-Alignment

Actuators & Adjusters

## Motion Control





### Micrometer Heads 1/4" (6.5mm) Travel



148-206

These small (0.25" range) micrometers are ideal for optical mounts. The spherical spindle tip provides a single point of contact, thus eliminating alignment errors. Thorlabs can integrate these micrometers into many of our standard products; call for details.

Range: 0 - 0.250" / .001" per Div.  
0 - 6.50mm / 10µm per Div.



ITEM#	\$	£	€	¥	RANGE
148-206	\$ 66.00	£ 49.50	€ 69,30	¥ 11,220	0-0.250"
148-205	\$ 66.00	£ 49.50	€ 69,30	¥ 11,220	0-6.5mm



### Micrometer Heads 1/2" (13mm) Travel

148-811



Range: 0 - 0.500" / .001" per Div.  
0 - 13.00mm / 10µm per Div.

These micrometer heads feature 1/2" (13mm) of travel and spherical spindle tips.




ITEM#	\$	£	€	¥	RANGE
148-811	\$ 67.00	£ 50.25	€ 70,35	¥ 11,390	0-0.500"
148-801	\$ 67.00	£ 50.25	€ 70,35	¥ 11,390	0-13mm

### Micrometer Heads 1.0" (25mm) Travel

150-811



Range: 0 - 1.000" / .001" per Div.  
0 - 25mm / 10µm per Div.



These micrometer heads feature 1" (25mm) of travel and spherical spindle tips.


ITEM#	\$	£	€	¥	RANGE
150-811	\$ 76.00	£ 57.00	€ 79,80	¥ 12,920	0-1.000"
150-801	\$ 76.00	£ 57.00	€ 79,80	¥ 12,920	0-25mm

### Micrometer Heads 2.0" (50mm) Travel

151-412S



Range: 0 - 2.000" / .001" per Div.  
0 - 50mm / 10µm per Div.



These micrometer heads feature 2" (50mm) of travel and spherical spindle tips.

ITEM#	\$	£	€	¥	RANGE
151-412S	\$ 197.00	£ 147.75	€ 206,85	¥ 33,490	0-2.000"
151-411S	\$ 197.00	£ 147.75	€ 206,85	¥ 33,490	0-50mm