

Have you
seen our...

Science
Desk
Accessories

See
page 35

Ø1.5" ScienceDesk Post Mount Adapter



Two TA2 Post Mount Adapters shown with PSY211 Overhead Shelf (Posts and Shelf Sold Separately)

- Mount a Ø1.5" ScienceDesk Post (PSY161) to a Breadboard or Optical Table
- Metric and Imperial Compatible
- For Mounting ScienceDesk Accessories to an Optical Table or Breadboard

The TA2 Post Adapter allows a Ø1.5" ScienceDesk post to be fitted to an optical table or breadboard. With this adapter, post-mounted ScienceDesk accessories such as those shown below can be used on almost any optical breadboard or table. The TA2 comes with 1/4" - 20 and M6 capscrews for securing it to the table surface. The ScienceDesk Post attaches to the adapter using the included M16 capscrew.



TA2

Post Mount Adapter

ITEM #	\$	£	€	RMB	DESCRIPTION
TA2	\$ 30.00	£ 21.60	€ 26,10	¥ 239.10	Ø1.5" ScienceDesk Post Mount Adapter
PSY161	\$ 95.00	£ 68.40	€ 82,65	¥ 757.15	ScienceDesk Post, Ø1.5" x 29.5" (Ø38 mm x 750 mm)

ScienceDesk Supported by Ø1.5" ScienceDesk Posts



PSY191

ScienceDesk Shelves are constructed using Melamine-faced MDF. The shelves presented here are designed to be supported with Ø1.5" ScienceDesk Posts (PSY161).




PSY211

Posts

ITEM #	\$	£	€	RMB	L X W	LOAD CAPACITY	PSY161 POSTS NEEDED
PSY211	\$ 330.00	£ 237.60	€ 287,10	¥ 2,630.10	300 mm x 900 mm	50 kg (110 lbs)	2 (Included)
PSY201/S*	\$ 160.00	£ 115.20	€ 139,20	¥ 1,275.20	300 mm x 900 mm	50 kg (110 lbs)	2 (Not Included)
PSY212	\$ 335.00	£ 241.20	€ 291,45	¥ 2,669.95	300 mm x 1200 mm	50 kg (110 lbs)	2 (Included)
PSY202/S*	\$ 170.00	£ 122.40	€ 147,90	¥ 1,354.90	300 mm x 1200 mm	50 kg (110 lbs)	2 (Not Included)
PSY191	\$ 177.00	£ 127.44	€ 153,99	¥ 1,410.69	300 mm x 278 mm	20 kg (40 lbs)	1 (Included)

*The /S option allows for the purchase of the shelf only (no posts).

ScienceDesk Computer Accessories Supported by Ø1.5" ScienceDesk Posts



PSY121

The monitor brackets are compatible with 75 mm x 75 mm and 100 mm x 100 mm VESA patterns on the back of most monitors. A Ø1.5 ScienceDesk Post is not included.

The keyboard holder comes with a Ø1.5 post just tall enough for the keyboard holder. The included post can be used with the TA2 ScienceDesk Post Adapter.



PSY132

Computer Accessories

ITEM #	\$	£	€	RMB	DESCRIPTION
PSY120	\$ 105.00	£ 75.60	€ 91,35	¥ 836.85	Flat Screen Support Bracket, 10 kg Load Capacity
PSY121	\$ 150.00	£ 108.00	€ 130,50	¥ 1,195.50	Flat Screen Support Bracket with Articulated Arm, 10 kg Load Capacity
PSY132	\$ 399.00	£ 287.28	€ 347,13	¥ 3,180.03	Keyboard Holder with Short Ø1.5" Mounting Post