

## Storage Boxes

Thorlabs now offers á la carte the same storage cases that protect the lenses and filters purchased in our kits. These storage boxes are perfect for any laboratory looking to organize and protect their optics investment.

**NEW**  
products



KT01

### Storage Box for Mounted Ø1" Optics (10 Slots: 0.45" Thick)

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT01	10	Ø1"	0.45" (11.4 mm)	Mounted in SM1 Lens Tubes	FGK01: Colored Glass Filter Kit (See Page 347) NEK01: Absorptive ND Filter Kit (See Page 348)



KT02

### Storage Box for Mounted Ø1/2" Optics

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT02	28	Ø1/2"	0.83" (21 mm)	Mounted in SM05 Lens Tubes	NEK02: Absorptive ND Filter Kit (See Page 348)



KT03

### Storage Box for Unmounted 2" Square Optics

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT03	10	2" x 2"	0.25" (6.4 mm)	Unmounted	FGK01S: Colored Glass Filter Kit (See Page 347) NEK01S: Absorptive ND Filter Kit (See Page 348)



KT04

### Storage Box for Mounted Ø1" Optics (26 Slots: 0.45" Thick, 9 Slots: 0.65" Thick)

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT04	35	Ø1"	9 Slots: 0.65" (16.5 mm) 26 Slots: 0.45" (11.4 mm)	Mounted in SM1 Lens Tubes	LSB01 Series: N-BK7 Plano-Convex Lens Kits (See Page 345)
					LSB02 Series: N-BK7 Bi-Convex Lens Kits (See Page 345)
					LSB03 Series: N-BK7 Plano- and Bi-Convex Lens Kits (See Page 345)
					LSB04 Series: N-BK7 Plano- and Bi-Convex/Concave Lens Kits (See Page 345)
					LSB05 Series: UV Fused Silica-Plano Convex Lens Kits (See Page 345)



KT05

### Storage Box for Mounted Ø1" Optics (35 Slots: 0.65" Thick)

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT05	35	Ø1"	0.65" (16.5 mm)	Mounted in SM1 Lens Tubes	LSB08 Series: Ø1" Achromatic Lens Kits (See Page 342)



KT06

### Storage Box for Mounted Ø2" Optics

ITEM#	CAPACITY	OPTIC SIZE	MAX THICKNESS	OPTIC HOUSING	OPTICS KITS USING THIS STORAGE BOX
KT06	10	Ø2"	Ø0.5" (12.7 mm)	Mounted in SM2 Lens Tubes	LSC01 Series: Ø2" Plano-Convex Lens Kits (See Page 345)

ITEM#	\$	£	€	RMB	DESCRIPTION
KT01	\$ 86.70	£ 60.20	€ 77.00	¥ 732.10	Storage Box for Mounted Ø1" Round Optics (Max. Capacity: 10)
KT02	\$ 86.70	£ 60.20	€ 77.00	¥ 732.10	Storage Box for Mounted Ø1/2" Round Optics (Max. Capacity: 28)
KT03	\$ 81.60	£ 56.60	€ 72.50	¥ 689.10	Storage Box for Unmounted 2" Square Optics (Max. Capacity: 10)
KT04	\$ 127.50	£ 88.40	€ 113.20	¥ 1,076.70	Storage Box for Mounted Ø1" Round Optics (26 Slots: 0.45" Thick; 9 Slots: 0.65" Thick)
KT05	\$ 127.50	£ 88.40	€ 113.20	¥ 1,076.70	Storage Box for Mounted Ø1" Round Optics (35 Slots: 0.65" Thick)
KT06	\$ 86.70	£ 60.20	€ 77.00	¥ 732.10	Storage Box for Mounted Ø2" Round Optics (Max. Capacity: 10)