

2" (50.8 mm) Travel, Differential Adjuster

€ 173,20

£ 135.20

This differential drive has been designed for use on single-axis stages equipped with a $\emptyset 1/2$ " ($\emptyset 12.7 \text{ mm}$) actuator mounting hole. The large-diameter differential drive actuator provides fine, ultra-smooth control over motion. The removable brass bushing allows this actuator to be used in any OEM application requiring 1" (25 mm) of overall travel.

0.55" - 1.06"

(14 - 27 mm)

RMB

¥ 1,646.60

1.46'

(37 mm)

DESCRIPTION

Differential Adjuster, 8 mm of Travel

Please refer to our website for complete models and drawings.



THORLABS

ITEM#

DM10A

Output via Non-Rotating Hard Steel Ball

\$

\$ 195.00